


PatientKeeper® Mobile Clinical Results™ on the iPhone®

If you have connectivity to your system's network, you can use your iPhone to review your patients' clinical results. To start, launch a browser from your iPhone, go to the URL listed below, then log in to PatientKeeper. Your patient list is displayed.


Select a patient from your patient list ... then select a Module (like Medications).



iPhone Browser URL: _____
 Support: _____


How much Detail do you want to see? Sometimes, a Summary is all you need, like a list of all the patient's current medications. Other times, you want to know all the details, like when a drug was ordered, or whether a dose was given.

Tap an item on the Summary ...



All of Molly Darr's recent medications are displayed. Let's look more closely at the Albuterol order.

... to see the Details.



The details for Albuterol show that the order was started on November 20.

Scroll down to the MAR section to see whether a dose was Given or Not Given.


Scroll down to see MAR data (if available).



How much Patient Information do you want to display on your device? The last 7 days of clinical data might be just the right amount of information for you. Or maybe you only need the past five days. You can leave the default settings in place, or change the settings any time.



You can control how much patient data you want to display on your device. From the go menu, tap Settings to display the Settings screens.




You can change the timeframe and module settings to show more or fewer results. Tap the dropdown list and choose a value. Tap Save to save the new values.

PatientKeeper® Mobile Clinical Results™

Use these Controls to Manage What's On Your Patient List


Manually add or remove a patient. Your patient list is always updated with your current and recent patients. You don't need to manually add or remove patients unless you want.



Sort the patient list. Sort alphabetically (A-Z or Z-A) by Name or by Location.





Select a Filter. Show All the patients on your list, or reduce the number of patients you see on your list with a filter.



Select a Profile. If your organization offers Profiles, you can choose a different Profile. Tap the go button and choose a profile. Your device syncs, and patients and data that match the profile are downloaded to your device.

Use a Graph to See How Results are Trending Over Time. When you review the details of a Lab or Vital, you can switch between a table view and a graphical view of the data.



The graph shows results trending over time. If you want to see the current value, normal range, and any comment, tap the table icon..

Value:	5.8 HH
Normal Range:	3.6-5.0 mEq/L
Status:	F
Comment:	Specimen not hemolyzed

Use the BACK button to move back a step



Vitals Summary. Tapping Back from the Summary list takes you back to your Patient List.

If you are on a Details screen, like this Vitals FSBG graph, Back takes you back to the Summary list.

Patient List. You always return to the Patient List when you click the Back button from a Summary page.